

DENOVO

Your Source for Quality Dental Products

(800) 854-7949

Or (626) 939-5000

FAX: (800) 847-8599 or (626) 939-5020

Internet: www.denovodental.com

email: sales@denovodental.com

General Purpose Molar Bands

Upper (Maxillary)	Size	Lower (Mandibular)
146	23	246
148	24	248
150	25	250
152	26	252
154	27	254
156	28	256
157	28½	257
158	29	258
159	29½	259
160	30	260
161	30½	261
162	31	262
163	31½	263
164	32	264
165	32½	265
166	33	266
167	33½	267
168	34	268
169	34½	269
170	35	270
171	35½	271
172	36	272
173	36½	273
174	37	274
175	37½	275
176	38	276
177	38½	277
178	39	278
179	39½	279
180	40	280
181	40½	281
182	41	282
183	41½	283
184	42	284

DENOVO general purpose molar bands are made in regular temper and narrow occlusal-gingival height. These bands, designed with maxillary and mandibular shapes, have universal anatomy to fit both 1st primary, 2nd primary, and 1st permanent molars. Each band is fully proportioned to provide a gradual increase in height as band size increases.

Sizes are available from 23-42, with ½ sizes available for precise fitting. The bands are completely festooned.

To use the order form at left:

1. Enter the quantities for each item in the white boxes (upper bands on the left, lower bands on right).
2. If you require complete part numbers for a purchase order, add the prefix 160 to all part numbers (an upper size 25 band has a part # of 160-150).
3. Quantity discounts are available for your total order. Add up all of the bands ordered to find out

Total Bands ordered	Price per band
1-29 bands	Retail
30-49 bands	5% discount
50-99 bands	10% discount
100 or more	15% discount

Molar Band Kits

Molar Band Kit Upper & Lower	161-100	
Molar Band Kit Maxillary (Upper)	161-101	
Molar Band Kit Mandibular (Lower)	161-102	

Kit Description	Part #
Maxillary & Mandibular Band Kit (Contains 288 Bands)	161-100
Maxillary Band Kit (Contains 144 Bands)	161-101
Mandibular Band Kit (Contains 144 Bands)	161-102

Name _____

Address _____

City, State, Zip _____

Phone _____

Comment _____