

DENOVO

Your Source for Quality Dental Products

(800) 854-7949

Or (626) 939-5000

FAX: (800) 847-8599 or (626) 939-5020

Internet: www.denovodental.com

email: sales@denovodental.com

Stainless Steel Crown Order Form

Stainless Steel Crowns

Size	1st Primary		2nd Primary		1st Permanent	
UR1	121		221		321	
UR2	122		222		322	
UR3	123		223		323	
UR4	124		224		324	
UR5	125		225		325	
UR6	126		226		326	
UR7	127		227		327	
UR8	n/a		n/a		328	
UL1	111		211		311	
UL2	112		212		312	
UL3	113		213		313	
UL4	114		214		314	
UL5	115		215		315	
UL6	116		216		316	
UL7	117		217		317	
UL8	n/a		n/a		318	
LR1	141		241		341	
LR2	142		242		342	
LR3	143		243		343	
LR4	144		244		344	
LR5	145		245		345	
LR6	146		246		346	
LR7	147		247		347	
LR8	n/a		n/a		348	
LL1	131		231		331	
LL2	132		232		332	
LL3	133		233		333	
LL4	134		234		334	
LL5	135		235		335	
LL6	136		236		336	
LL7	137		237		337	
LL8	n/a		n/a		338	

Crowns are pre-trimmed for proper tooth coverage; anatomical gingival contour saves chair time when fitting crown.

Parallel crown walls saves chair time, minimizing the need for buccal, lingual, and mesial tooth reduction.

Consistent wall thickness helps prevent bite through and minimizes trial fitting distortion.

Shallow occlusal anatomy reduces rocking during mastication, with minimal occlusal interference.

DENOVO Stainless Steel Molar Crowns allow you to save healthy tooth structure, and eliminate time-consuming and unnecessary over-reduction of tooth surface. Shallow occlusal surface provides confidence of achieving optimum functional occlusion.

Stainless Steel Crown Kits

Primary kits contain 56 crowns.

The 1st Permanent kit contains 64 crowns.

Crown Kit	Part #
1st Primary Kit	902-100
2nd Primary Kit	902-200
1st Permanent Kit	902-300

Crown Quantity Discounts

Total Order	Price
1-24 Crowns	Retail
25-49 Crowns	5% off
50-99 Crowns	10% off
100+ Crowns	15% off

Name _____

Address _____

City, State, Zip _____

Phone _____

Comment _____